



PC-80A

• Easy ECG Monitor

The Scope of Application:

PC-80A Easy ECG Monitor is designed for home use with the purpose of health care, and is suitable for users who may suffer from cardiovascular disease, such as hypertension, diabetes, coronary heart disease, etc., or those who have the symptoms of heart-throb, stuffy chest. It can be used to check their heart conditions regularly, and the stored ECG records can be used as a reference for clinical ECG examination.

Features:

- Small, portable and easy to operate
- 24 records storage and review, each record with 30-second ECG waveform and measurement result
- Measurement of one channel ECG anytime, anywhere you like
- Display of ECG waveform, heart rate, measurement results and battery status on the backlight LCD
- Auto power-off while no key is pressed 35 seconds
- Lock function to data records to avoid inadvertent deletion
- 2 batteries can take 500 times measurement



CE 1011

Your health, our concern!