



PC-68

Wrist Oximeter •



The Scope of Application:

The PC-68 Wrist Oximeter is small, lightweight, and available for home use as well as clinics in long-term and comfortable monitoring. It can be simply worn on the patient's wrist as the watch to start the monitoring conveniently. This device is applicable to overnight monitoring with sleep or respiration monitoring devices, and can be used in sports to measure SpO₂ as well as pulse. The advanced technology ensures a reliable and accurate measurement.

Features:

- Automatic power-on function
- Color OLED to display SpO₂, pulse rate, and P1 as well as the bar-graph of pulse
- Stored measurement results transmission to PC by serial port
- Two AAA batteries are required for operation with low power consumption
- Battery voltage indication
- Alarm function, high & low alarm limits adjustable
- Time & date display, adjustable by the user
- Time setting for the automatic measurement
- Storage mode selectable: 1 group/s, 1 group/4s, and 1 group/8s



Your health, our passion!

Note: Specifications subject to change without prior notice